



NYLON

NYLON

Non-absorbable, synthetic surgical monofilament thread




- Polyamide
- Uncoated
- Thread may be blue or black or undyed
- Size 10/0 – 5 USP
- Features
 - Uniform structure with smooth surface
 - Flexible and elastic tightening of knots
 - Very good bindings properties
 - High durability of knot
 - Tensile strength
- Indications
 - Soft tissues (intra-dermal and dermal suturing)
 - General, Pediatric, Plastic, Vascular surgery
 - Urology
 - Ophthalmology
- Period of validity 5 years

yavo

YAVO Sp. z o. o. 97-400 Bełchatów, ul. Bawełniana 17, Polska
tel.: +48 44 632 94 11, fax:+48 44 632 49 17,
e-mail: info@yavo.com.pl, www.yavo.com.pl



NYLON



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament




spatula (lanced needle)

| length | needles | | thread length (cm) | 10/0 | 9/0 | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 | |
|-------------------|-----------|----------|--------------------|----------|--------|--------|--------|--------|-----|-----|--|
| | curvature | diameter | | | | | | | | | |
| 6 mm | 3/8 | 150 µm | | 7 black | LZ63PA | | | | | | |
| | | | | 15 black | LZ55PA | | | | | | |
| | | | | 20 black | LZ37PA | | | | | | |
| 6 mm micropoint | 3/8 | | | 30 black | LZ60PA | | | | | | |
| 6 mm micropoint | 3/8 | 150 µm | | 30 black | LZ52PA | LZ54PA | | | | | |
| 6 mm | 3/8 | 150 µm | | 30 black | LZ10PA | LZ11PA | | | | | |
| 6 mm | 3/8 | 200 µm | | 30 black | LZ12PA | LZ13PA | LZ26PA | LZ76PA | | | |
| | | | | 45 black | | | LZ57PA | | | | |
| 6 mm | 3/8 | 200 µm | | 30 black | LZ14PA | LZ15PA | | | | | |
| | | | | 45 black | LZ61PA | | | | | | |
| 6,19 mm | 3/8 | 150 µm | | 30 black | LZ49PA | LZ56PA | | | | | |
| 6,19 mm | 3/8 | 200 µm | | 30 black | LZ18PA | LZ19PA | | | | | |
| 6,19 mm | 3/8 | 200 µm | | 30 black | LZ20PA | LZ21PA | | | | | |
| 6,2 mm | 3/8 | | | 30 black | LZ30PA | LZ31PA | | | | | |
| 6,2 mm micropoint | 3/8 | 200 µm | | 30 black | LZ78PA | | | | | | |
| 6,4 mm | 3/8 | | | 30 black | | | LZ41PA | | | | |
| 6,4 mm | 3/8 | 200 µm | | 30 black | LZ22PA | LZ23PA | LZ24PA | | | | |
| 6,5 mm | 3/8 | 150 µm | | 30 black | | LZ67PA | | | | | |
| 6,5 mm | 3/8 | 220 µm | | 30 black | LZ35PA | LZ36PA | | | | | |
| 6,5 mm micropoint | 3/8 | 220 µm | | 30 black | LZ64PA | LZ65PA | | | | | |
| | | | | 45 black | LZ79PA | | | | | | |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



spatula (lanced needle)

| length | needles | | curvature | diameter | thread length (cm) | color | 10/0 | 9/0 | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 |
|----------------------|-----------|----------|-----------|----------|--------------------|-------|--------|--------|--------|-----|--------|--------|-----|
| | curvature | diameter | | | | | | | | | | | |
| 6,5 mm micropoint | 3/8 | 150 µm | | | 30 | black | LZ69PA | | | | | | |
| | | | | | 45 | black | | | LZ77PA | | | | |
| 8 mm | 3/8 | | | | 30 | black | | | | | LZ40PA | | |
| 5,5 mm | 1/2 | | | | 30 | black | LZ28PA | | | | | | |
| 5,5 mm micropoint | 1/2 | 150 µm | | | 30 | black | LZ53PA | | | | | | |
| 5,6 mm | 1/2 | 150 µm | | | 20 | black | LZ03PA | LZ04PA | | | | | |
| 5,6 mm | 1/2 | 150 µm | | | 30 | black | LZ06PA | LZ07PA | LZ08PA | | | | |
| 6 mm | 1/2 | 150 µm | | | 15 | black | LZ33PA | | | | | | |
| | | | | | 30 | black | LZ16PA | | LZ59PA | | | | |
| 6 mm | 1/2 | 150 µm | | | 30 | black | LZ17PA | | LZ39PA | | | | |
| 6 mm micropoint | 1/2 | | | | 30 | black | | | LZ70PA | | | | |
| 6,5 mm micropoint | 1/2 | 220 µm | | | 30 | black | LZ74PA | | | | | | |
| 8 mm | 1/4 | 280 µm | | | 45 | blue | | | | | | LZ75PA | |

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

| length | needles curvature | thread length (cm) color | 10/0 | 9/0 | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 |
|--------|----------------------|--------------------------------|---------|---------|---------|---------|---------|---------|---------|
| | | | | | | | | | |
| 4 mm | 3/8 | 13 black | RS115PA | RS248PA | | | | | |
| | | 15 black | RS77PA | | | | | | |
| 5 mm | 3/8 | 13 black | RS212PA | RS214PA | RS224PA | | | | |
| | | 15 black | RS130PA | | | | | | |
| | | 45 black | RS200PA | | | | | | |
| 6,4 mm | 3/8 | 13 black | | | RS245PA | | | | |
| | | 15 black | | RS02PA | RS03PA | | | | |
| 6,4 mm | 3/8 | 45 black | RS89PA | RS04PA | RS05PA | | | | |
| 8 mm | 3/8 | 60 blue | | | | RS182PA | | | |
| 9,5 mm | 3/8 | 75 blue | | | | RS08PA | RS09PA | | |
| 10 mm | 3/8 | 45 blue | | | RS201PA | | RS10PA | | |
| 11 mm | 3/8 | 60 blue | | | | RS106PA | | | |
| 12 mm | 3/8 | 45 blue | | | | | RS11PA | RS12PA | |
| 12 mm | 3/8 | 45 blue | | | | | RS79PA | | |
| 13 mm | 3/8 | 60 blue | | | | | RS13PA | | |
| | | 75 blue | | | | RS158PA | RS159PA | RS184PA | RS185PA |
| 18 mm | 3/8 | 90 blue | | | | | | RS160PA | |

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

| length | needles curvature | thread length (cm) color | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 |
|--------|----------------------|--------------------------------|-------|-----|------------|---------|---------|---------|---------|
| | | | 19 mm | 3/8 | 45 blue | | | | RS156PA |
| 22 mm | 3/8 | 45 blue | | | | RS194PA | | | |
| 22 mm | 3/8 | 90 blue | | | | | RS161PA | | |
| 24 mm | 3/8 | 45 blue | | | | | RS204PA | RS121PA | |
| | | 75 blue | | | | | RS125PA | RS124PA | RS123PA |

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

| length | needles curvature | thread length (cm) color | 10/0 | 9/0 | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 |
|--------|----------------------|--------------------------------|----------------|-----|-----|---------------|---------------|----------------|-----|
| | | | | | | | | | |
| 5 mm | 1/2 | 13 black | RS231PA | | | | | | |
| 10 mm | 1/2 | 75 blue | | | | RS85PA | | | |
| 12 mm | 1/2 | 75 blue | | | | | RS83PA | RS84PA | |
| 13 mm | 1/2 | 75 blue | | | | | | RS218PA | |

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

| length | needles curvature | | thread length (cm) color | 6/0 5/0 4/0 3/0 2/0 0 1 | | | | | | |
|----------------|----------------------|--|--------------------------------|-------------------------|---------|---------|---------|---------|---------|---------|
| | | | | | | | | | | |
| 16 mm | 1/2 | | 75 blue | | RS230PA | RS82PA | RS227PA | | | |
| 16 mm | 1/2 | | 75 blue | RS171PA | | RS132PA | | | | |
| | | | 90 blue | | RS14PA | RS15PA | RS205PA | | | |
| 17 mm | 1/2 | | 75 blue | | | | RS193PA | | | |
| 20 mm | 1/2 | | 75 blue | | RS16PA | RS17PA | RS18PA | RS19PA | RS20PA | |
| 20 mm | 1/2 | | 90 blue | | | RS105PA | | | | |
| | | | | | | | | | | |
| 22 mm | 1/2 | | 75 blue | | | RS236PA | RS237PA | | | |
| | | | 90 blue | | | RS167PA | RS180PA | RS150PA | | RS149PA |
| 26 mm | 1/2 | | 75 blue | | | | RS24PA | RS25PA | RS97PA | RS127PA |
| | | | 90 blue | | | | | | RS26PA | |
| 26 mm | 1/2 | | 90 blue | | | | | RS113PA | | |
| | | | 120 blue | | | | | RS104PA | | |
| 30 mm | 1/2 | | 75 blue | | RS107PA | RS155PA | | RS32PA | RS33PA | RS34PA |
| | | | 90 blue | | | | | RS35PA | RS142PA | RS143PA |
| | | | 100 blue | | | | | RS36PA | RS37PA | RS38PA |
| 30 mm heavy | 1/2 | | 75 blue | | | | | RS102PA | | |
| 30 mm | 1/2 | | 150 blue | | | | | RS223PA | RS153PA | |
| 35 mm | 1/2 | | 150 blue | | | | | | | RS137PA |
| 36 mm heavy | 1/2 | | 45 blue | | | | | | RS170PA | RS169PA |

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament




round needle

| length | needles curvature | thread length (cm) color | thread | | | | | | | |
|----------------|----------------------|--------------------------------|--------|-----|-----|--------|--------|---------|---------|---------|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
| 37 mm | 1/2 | 75 blue | | | | RS53PA | RS54PA | RS55PA | RS128PA | |
| | | 90 blue | | | | | | | RS56PA | |
| | | 100 blue | | | | | | RS57PA | RS76PA | |
| 40 mm | 1/2 | 75 blue | | | | | | RS60PA | RS61PA | RS62PA |
| | | 90 blue | | | | | | | | RS110PA |
| | | 100 blue | | | | | | RS63PA | RS64PA | RS65PA |
| 40 mm heavy | 1/2 | 75 blue | | | | | | | | RS221PA |
| | | 100 blue | | | | | | RS202PA | | |
| 40 mm | 1/2 | 150 blue | | | | | | RS139PA | | |
| | | 200 blue | | | | | | RS238PA | | |
| 40 mm heavy | 1/2 | 150 blue | | | | | | RS203PA | | |
| 48 mm | 1/2 | 90 blue | | | | | | | | RS108PA |
| | | 100 blue | | | | | | RS189PA | RS66PA | RS67PA |
| | | 200 blue | | | | | | | RS165PA | |
| | | 240 blue | | | | | | | RS186PA | |
| 48 mm heavy | 1/2 | 100 blue | | | | | | RS247PA | | |
| 48 mm | 1/2 | 150 blue | | | | | | RS190PA | RS181PA | |
| | | 240 blue | | | | | | | RS188PA | |
| 48 mm heavy | 1/2 | 150 blue | | | | | | RS88PA | | RS103PA |
| | | 200 blue | | | | | | | RS220PA | |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament




round needle

| length | needles curvature | thread length (cm) color | thread | | | | | | | 1 | 2 |
|----------------|----------------------|--------------------------------------|--------|-----|-----|-----|---|--------|---------|---------|---------|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | | |
| 50 mm | 1/2 | 100 blue | | | | | | | | RS68PA | RS234PA |
| | | 150 blue | | | | | | RS69PA | RS70PA | | |
| 50 mm heavy | 1/2 | 100 blue | | | | | | | | RS196PA | |
| 50 mm | 1/2 | 150 blue | | | | | | | | RS233PA | |
| | | 240 blue | | | | | | | | RS114PA | |
| 60 mm | 1/2 | 75 blue 2 vinyl tubes 60 mm | | | | | | | RS242PA | | |
| | | 90 blue | | | | | | | | RS241PA | |
| 65 mm | 1/2 | 150 blue | | | | | | | | RS71PA | RS72PA |
| 65 mm | 1/2 | 120 blue | | | | | | | RS222PA | | |
| | | 150 blue | | | | | | | | | RS192PA |
| | | 240 blue | | | | | | | RS187PA | | |
| 76 mm | 1/2 | 75 blue | | | | | | | RS100PA | RS73PA | |
| 80 mm | 1/2 | 100 blue | | | | | | | | RS166PA | |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

| length | needles curvature | thread length (cm) color | | | | | | | | |
|--------|----------------------|--------------------------------|-----|-----|-----|---------|---------|---------|---|--|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
| 60 mm | straight | 75 blue | | | | RS148PA | RS91PA | | | |
| | | 90 blue | | | | | RS129PA | RS154PA | | |
| 60 mm | straight | 90 blue | | | | | RS249PA | | | |
| 65 mm | straight | 90 blue | | | | | RS229PA | | | |
| 70 mm | straight | 75 blue | | | | | RS199PA | | | |
| 70 mm | straight | 75 blue | | | | | RS98PA | | | |

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

| length | needles curvature | | thread length (cm) color | 6/0 5/0 4/0 3/0 2/0 0 1 | | | | | | | |
|--------|----------------------|--|--------------------------------|-------------------------|---------|---------|---------|---------|---------|---------|--|
| | | | | | | | | | | | |
| 8 mm | 3/8 | | 45 blue | DK130PA | | | | | | | |
| 12 mm | 3/8 | | 45 blue | DK71PA | DK95PA | DK67PA | DK85PA | DK86PA | | | |
| | | | 75 blue | DK170PA | DK127PA | DK68PA | | | | | |
| 12 mm | 3/8 | | 75 blue | | DK173PA | | | | | | |
| 16 mm | 3/8 | | 45 blue | DK01PA | DK02PA | DK03PA | DK87PA | DK88PA | | | |
| | | | 45 black | DK174PA | | | | | | | |
| | | | 75 blue | DK04PA | DK05PA | DK06PA | DK128PA | | | | |
| 16 mm | 3/8 | | 45 blue | DK07PA | DK08PA | DK168PA | | | | | |
| | | | 75 blue | DK09PA | DK10PA | | | | | | |
| 19 mm | 3/8 | | 45 blue | DK123PA | DK11PA | DK12PA | DK13PA | DK89PA | | | |
| | | | 75 blue | | DK14PA | DK15PA | DK16PA | DK129PA | | | |
| | | | 90 blue | | DK17PA | DK18PA | DK19PA | DK20PA | | | |
| 19 mm | 3/8 | | 45 blue | | DK132PA | DK107PA | DK92PA | | | | |
| 24 mm | 3/8 | | 45 blue | | DK100PA | DK21PA | DK22PA | DK23PA | | | |
| | | | 75 blue | | | DK24PA | DK25PA | DK26PA | DK79PA | | |
| | | | 90 blue | | | DK27PA | DK28PA | DK29PA | DK172PA | | |
| 24 mm | 3/8 | | 75 blue | | | DK154PA | DK146PA | | | | |
| 26 mm | 3/8 | | 45 blue | | | DK30PA | DK31PA | DK32PA | | | |
| | | | 75 blue | | | DK33PA | DK34PA | DK35PA | DK69PA | DK165PA | |
| | | | 90 blue | | | DK36PA | DK37PA | DK38PA | DK115PA | | |
| 26 mm | 3/8 | | 45 blue | | | | DK156PA | | | | |
| | | | 75 blue | | | | DK114PA | DK112PA | | | |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.












NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

| length | needles curvature | thread length (cm) color | | | | | | | |
|-------------------|---|--------------------------------|-----|--------|---------|---------|---------|---------|---------|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 |
| 30 mm | 3/8  | 45 blue | | | DK135PA | DK39PA | DK40PA | | |
| | | 75 blue | | | DK41PA | DK42PA | DK43PA | DK75PA | |
| | | 90 blue | | | DK72PA | DK73PA | DK74PA | DK145PA | |
| 35 mm | 3/8  | 75 blue | | | DK133PA | DK103PA | | | |
| | | 90 blue | | | DK163PA | DK84PA | DK105PA | | |
| 36 mm cosmetic | 3/8  | 75 blue | | | DK113PA | | | | |
| 39 mm | 3/8  | 45 blue | | | | | | DK158PA | |
| | | 75 blue | | DK93PA | DK94PA | DK78PA | DK46PA | DK47PA | |
| | | 90 blue | | | | DK122PA | DK137PA | DK142PA | DK157PA |
| 40 mm | 3/8  | 45 blue | | | | DK48PA | DK49PA | | |
| | | 75 blue | | | | DK50PA | DK51PA | DK101PA | |
| | | 90 blue | | | | DK150PA | | | |
| 40 mm heavy | 3/8  | 45 blue | | | | DK121PA | | | |
| | | 75 blue | | | | DK125PA | | | |
| 48 mm | 3/8  | 100 blue | | | | DK140PA | | | |
| 60 mm | 3/8  | 100 blue | | | | DK99PA | DK76PA | | DK77PA |
| 90 mm | 3/8  | 45 blue | | | | | | DK124PA | |
| | | 75 blue | | | | | DK153PA | | |
| | | 90 blue | | | | | | DK59PA | DK60PA |
| | | 100 blue | | | | DK61PA | DK62PA | DK63PA | DK64PA |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

| length | needles curvature | thread length (cm) color | 5/0 4/0 3/0 2/0 0 1 2 | | | | | | | |
|--------|---------------------------|--------------------------------|-----------------------|-----|---------|---------|--------|---------|---------|---------|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
| 90 mm | 3/8 2 vinyl tubes 40mm | 2x50 blue | | | | | | | | DK155PA |
| 90 mm | 3/8 | 75 blue | | | | | | | | DK166PA |
| 20 mm | 1/2 | 45 blue | | | DK181PA | DK182PA | | | | |
| 39 mm | 1/2 | 75 blue | | | | DK136PA | DK97PA | | | |
| | | 90 blue | | | | | DK44PA | DK45PA | | |
| 40 mm | 1/2 | 45 blue | | | | | | DK111PA | | |
| | | 75 blue | | | | | | DK110PA | | |
| | | 90 blue | | | | | | | DK187PA | |
| | | 100 blue | | | | | | DK54PA | DK180PA | |
| | | 150 blue | | | | | | | DK188PA | |
| 48 mm | 1/2 | 75 blue | | | | | | | DK57PA | |
| | | 100 blue | | | | | | | DK58PA | |
| 90 mm | 1/2 | 100 blue | | | | | | | | DK65PA |

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

| length | needles curvature | thread length (cm) color | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
|--------|----------------------|--------------------------------|-----|-----|---------|---------|---------|---------|---|--|
| | | | | | | | | | | |
| 51 mm | straight | 75 blue | | | | DG030PA | DG031PA | DG032PA | | |
| | | 90 blue | | | | | DG040PA | | | |
| 55 mm | straight | 75 blue | | | | | DG07PA | | | |
| | | 100 blue | | | | | DG08PA | | | |
| 60 mm | straight | 75 blue | | | DG020PA | DG018PA | DG029PA | | | |
| | | 90 blue | | | | DG021PA | DG041PA | | | |
| 60 mm | straight | 90 blue | | | | | DG26PA | | | |
| 75 mm | straight | 75 blue | | | | DG019PA | DG043PA | | | |

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | | thread length (cm) color | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 | |
|-------------------|----------------------|--|--------------------------------|---------|----------|----------|----------|----------|----------|----------|----------|
| | | | | | | | | | | | |
| 8 mm | 3/8 | | 45 blue | DKO01PA | | | DKO226PA | DKO244PA | | | |
| 8 mm cosmetic | 3/8 | | 45 blue | | | | DKO227PA | | | | |
| | | | 45 black | | | | DKO379PA | DKO381PA | | | |
| 11 mm | 3/8 | | 45 blue | | DKO04PA | DKO05PA | | | | | |
| | | | 45 black | | | | DKO312PA | | | | |
| | | | 75 black | | | | | DKO206PA | | | |
| 11 mm cosmetic | 3/8 | | 45 blue | | DKO194PA | DKO235PA | DKO323PA | | | | |
| 12 mm | 3/8 | | 45 blue | | | | DKO08PA | DKO09PA | DKO10PA | DKO327PA | |
| | | | 75 blue | | | | DKO11PA | DKO12PA | DKO13PA | | |
| 12 mm cosmetic | 3/8 | | 45 blue | | | | DKO14PA | DKO15PA | DKO16PA | | |
| | | | 45 black | | | | DKO318PA | DKO317PA | | | |
| | | | 75 blue | | | | DKO17PA | DKO18PA | DKO19PA | | |
| 13 mm | 3/8 | | 45 blue | | | | DKO20PA | DKO21PA | DKO224PA | | |
| | | | 75 blue | | | | | DKO278PA | | | |
| 13 mm | 3/8 | | 45 blue | | | | | DKO346PA | | | |
| 13 mm cosmetic | 3/8 | | 45 blue | | | | DKO213PA | DKO22PA | DKO251PA | | |
| 16 mm | 3/8 | | 45 blue | | | | DKO23PA | DKO24PA | DKO25PA | DKO26PA | |
| | | | 45 black | | | | | DKO361PA | | | |
| 16 mm cosmetic | 3/8 | | 75 blue | | | | DKO27PA | DKO28PA | DKO29PA | DKO30PA | DKO220PA |
| | | | 45 blue | | | | DKO31PA | DKO32PA | DKO33PA | DKO333PA | |
| | | | 75 blue | | | | DKO34PA | DKO35PA | DKO36PA | | |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | thread length (cm) color | 6/0 5/0 4/0 3/0 2/0 0 1 | | | | | | |
|-------------------|----------------------|--------------------------------|-------------------------|----------|----------|----------|----------|----------|----------|
| | | | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 |
| 19 mm | 3/8 | 45 blue | DKO38PA | DKO39PA | DKO40PA | DKO41PA | DKO42PA | DKO358PA | |
| | | 45 black | DKO376PA | DKO360PA | | | | | |
| | | 75 blue | | DKO43PA | DKO44PA | DKO45PA | DKO46PA | DKO331PA | |
| | | 75 black | | DKO377PA | | | | | |
| | | 90 blue | | | DKO47PA | DKO48PA | DKO49PA | | |
| 19 mm cosmetic | 3/8 | 45 blue | | DKO50PA | DKO51PA | DKO52PA | DKO277PA | | |
| | | 45 black | | DKO316PA | | | | | |
| | | 75 blue | DKO53PA | DKO54PA | DKO55PA | DKO56PA | | | |
| 22 mm | 3/8 | 45 blue | | DKO279PA | DKO321PA | | | | |
| | | 75 blue | | | DKO276PA | DKO275PA | DKO263PA | | |
| | | 90 blue | | | | DKO192PA | | | |
| 24 mm | 3/8 | 45 blue | DKO234PA | DKO59PA | DKO60PA | DKO61PA | DKO62PA | DKO63PA | DKO301PA |
| | | 75 blue | | DKO64PA | DKO65PA | DKO66PA | DKO67PA | DKO68PA | |
| | | 90 blue | | | DKO69PA | DKO70PA | DKO71PA | | |
| 24 mm cosmetic | 3/8 | 45 blue | | DKO236PA | | DKO72PA | DKO73PA | | |
| | | 75 blue | | | | DKO74PA | DKO75PA | | |
| 26 mm | 3/8 | 45 blue | | DKO243PA | DKO189PA | DKO77PA | DKO78PA | DKO79PA | |
| | | 75 blue | | | DKO80PA | DKO81PA | DKO82PA | DKO83PA | DKO84PA |
| | | 90 blue | | | DKO285PA | | DKO320PA | DKO292PA | |
| | | 100 blue | | | | | | | DKO85PA |
| 26 mm cosmetic | 3/8 | 45 blue | | | DKO382PA | DKO201PA | DKO294PA | | |
| | | 75 blue | | | DKO282PA | DKO199PA | DKO200PA | | |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.












NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | thread length (cm) color | thread | | | | | | 1 | 2 |
|----------------|---|--------------------------------|----------|----------|----------|----------|----------|----------|----------|---|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | | | |
| 30 mm | 3/8  | 45 blue | DKO373PA | DKO86PA | DKO87PA | DKO88PA | DKO190PA | DKO313PA | | |
| | | 75 blue | | DKO90PA | DKO91PA | DKO92PA | DKO93PA | DKO94PA | DKO215PA | |
| | | 90 blue | | DKO256PA | DKO95PA | DKO96PA | DKO97PA | DKO98PA | | |
| 35 mm | 3/8  | 45 blue | | | DKO105PA | | | | | |
| | | 75 blue | | | DKO106PA | DKO107PA | DKO205PA | DKO265PA | | |
| | | 90 blue | | | DKO108PA | DKO109PA | DKO203PA | | | |
| 37 mm | 3/8  | 45 blue | | | DKO110PA | | | | | |
| | | 75 blue | | | DKO289PA | DKO219PA | DKO111PA | DKO112PA | | |
| | | 90 blue | | | DKO370PA | | | DKO113PA | | |
| 39 mm | 3/8  | 45 blue | | | DKO126PA | DKO208PA | | | | |
| | | 75 blue | | | DKO127PA | DKO128PA | DKO129PA | DKO130PA | DKO369PA | |
| | | 90 blue | | | DKO131PA | DKO132PA | DKO133PA | DKO134PA | DKO135PA | |
| | | 100 blue | | | | | DKO286PA | | | |
| 40 mm | 3/8  | 45 blue | | | | DKO204PA | DKO139PA | | | |
| | | 75 blue | | | DKO290PA | DKO140PA | DKO141PA | DKO142PA | DKO143PA | |
| | | 100 blue | | | | DKO144PA | DKO145PA | DKO271PA | DKO191PA | |
| 40 mm heavy | 3/8  | 75 blue | | | DKO384PA | DKO341PA | | DKO303PA | | |
| 40 mm cosmetic | 3/8  | 75 blue | | | | DKO221PA | | | | |
| 45 mm | 3/8  | 75 blue | | | | DKO157PA | | DKO268PA | | |
| | | 90 blue | | | | DKO158PA | DKO306PA | | | |
| | | 100 blue | | | | DKO160PA | DKO267PA | | | |
| 45 mm cosmetic | 3/8  | 75 blue | | | | | DKO385PA | | | |

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



DKO

NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | thread length (cm) color | | | | | | | | |
|--------|----------------------|--------------------------------|-----|-----|-----|-----|----------|----------|----------|----------|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
| 48 mm | 3/8 | 90 blue | | | | | DKO310PA | | DKO272PA | DKO240PA |
| | | 100 blue | | | | | DKO319PA | | DKO283PA | |
| 60 mm | 3/8 | 100 blue | | | | | DKO187PA | DKO171PA | DKO172PA | DKO173PA |
| 77 mm | 3/8 | 45 blue | | | | | | | | DKO174PA |
| | | 75 blue | | | | | | | DKO264PA | DKO175PA |
| | | 90 blue | | | | | | | DKO375PA | |
| 77 mm | 3/8 | 75 blue | | | | | | | DKO238PA | DKO176PA |
| 90 mm | 3/8 | 2x50 blue | | | | | | | DKO326PA | |
| | | 75 blue | | | | | | | | DKO269PA |
| | | 90 blue | | | | | | DKO178PA | | DKO179PA |
| | | 100 blue | | | | | DKO180PA | DKO181PA | DKO182PA | DKO183PA |

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | | thread length (cm) color | thread | | | | | | | | |
|------------------|----------------------|--|--------------------------------|----------|-----|-----|----------|----------|----------|----------|----------|--|
| | | | | 8/0 | 7/0 | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 | | |
| 8 mm | 1/2 | | 45 blue | DKO02PA | | | | DKO228PA | | | | |
| 8 mm cosmetic | 1/2 | | 45 blue | DKO295PA | | | | | | | | |
| 16 mm | 1/2 | | 45 blue | | | | DKO305PA | DKO366PA | DKO362PA | DKO367PA | | |
| 19 mm | 1/2 | | 45 blue | | | | | | DKO262PA | | | |
| | | | 90 blue | | | | | | DKO186PA | DKO239PA | | |
| 22 mm | 1/2 | | 75 blue | | | | | | | DKO354PA | DKO355PA | |
| | | | 90 blue | | | | | | | DKO58PA | | |
| 24 mm | 1/2 | | 75 blue | | | | | | | DKO237PA | DKO241PA | |
| | | | 90 blue | | | | | | | DKO76PA | | |
| 25 mm | 1/2 | | 75 blue | | | | | | DKO325PA | | | |

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | thread length (cm) color | | | | | | | | |
|----------------|----------------------|--------------------------------|-----|-----|-----|----------|----------|----------|----------|----------|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
| 26 mm | 1/2 | 75 blue | | | | DKO252PA | DKO254PA | | | DKO247PA |
| 30 mm | 1/2 | 75 blue | | | | DKO260PA | DKO99PA | DKO100PA | DKO101PA | DKO210PA |
| | | 90 blue | | | | | | DKO343PA | | |
| 37 mm | 1/2 | 75 blue | | | | | DKO115PA | DKO116PA | DKO117PA | DKO249PA |
| | | 90 blue | | | | | DKO118PA | DKO222PA | DKO119PA | |
| | | 100 blue | | | | | | | DKO120PA | |
| 40 mm | 1/2 | 75 blue | | | | | | DKO149PA | DKO150PA | DKO151PA |
| | | 90 blue | | | | | | DKO233PA | DKO152PA | DKO153PA |
| | | 100 blue | | | | | | DKO154PA | DKO155PA | DKO156PA |
| 40 mm heavy | 1/2 | 90 blue | | | | | DKO300PA | | | |
| 48 mm | 1/2 | 75 blue | | | | | | | DKO162PA | |
| | | 90 blue | | | | | | | DKO163PA | |
| | | 100 blue | | | | DKO364PA | DKO322PA | DKO164PA | DKO242PA | |
| 50 mm | 1/2 | 75 blue | | | | | | | DKO165PA | |
| | | 90 blue | | | | | | | DKO372PA | DKO332PA |
| | | 100 blue | | | | | | | DKO166PA | |
| 60 mm | 1/2 | 100 blue | | | | DKO253PA | DKO344PA | | | |

notes



DKO

NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



reverse cutting needle

| length | needles curvature | thread length (cm) color | thread | | | | | | | |
|--------|----------------------|--------------------------------|--------|---------|---------|---------|---------|---------|---|--|
| | | | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 | |
| 45 mm | straight | 75 blue | | | DG045PA | | | | | |
| 51 mm | straight | 90 blue | | DG02PA | DG03PA | DG04PA | DG05PA | DG06PA | | |
| 55 mm | straight | 100 blue | | | | DG09PA | | | | |
| 60 mm | straight | 75 blue | | DG033PA | DG010PA | DG011PA | | | | |
| | | 90 blue | | DG044PA | DG012PA | DG013PA | DG014PA | DG015PA | | |
| | | 100 blue | | | | DG016PA | DG017PA | | | |
| 60 mm | straight | 75 blue | | DG24PA | DG23PA | DG22PA | | | | |
| | | 90 blue | | | | | DG042PA | | | |
| 75 mm | straight | 75 blue | | | | DG039PA | DG028PA | | | |

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



**cosmetic reverse cutting needle
with a precise blade (micro point)**

| length | needles curvature | | thread length (cm) color | 6/0 5/0 4/0 3/0 2/0 0 1 | | | | | | | |
|--------|----------------------|--|--------------------------------|-------------------------|--------|--------|--------|--------|--------|---|--|
| | | | | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | |
| 13 mm | 3/8 | | 45 blue | ES23PA | ES13PA | ES29PA | | | | | |
| 16 mm | 3/8 | | 45 blue | ES18PA | ES05PA | ES06PA | ES24PA | | | | |
| | | | 75 blue | | ES19PA | ES25PA | | | | | |
| 19 mm | 3/8 | | 45 blue | | ES07PA | ES01PA | ES22PA | | | | |
| | | | 45 black | | ES04PA | | | | | | |
| | | | 75 blue | | | ES08PA | ES09PA | | | | |
| 24 mm | 3/8 | | 75 blue | | | ES21PA | ES10PA | ES11PA | | | |
| 26 mm | 3/8 | | 45 blue | | | | ES03PA | | | | |
| | | | 75 blue | | | | ES12PA | ES02PA | | | |
| 30 mm | 3/8 | | 75 blue | | | | ES32PA | ES14PA | ES27PA | | |
| 40 mm | 3/8 | | 75 blue | | | | | ES26PA | | | |

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament




round needle with a cutting edge (sharpened)

| length | needles curvature | thread length (cm) | color | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 |
|-------------|-------------------|--------------------|-------|--------|--------|--------|--------|--------|--------|--------|
| | | | | | | | | | | |
| 24 mm | 3/8 | 75 | blue | | DS16PA | DS17PA | DS18PA | | | |
| 26 mm | 3/8 | 75 | blue | | | | DS09PA | | | |
| 16 mm | 1/2 | 90 | blue | DS10PA | | | | | | |
| 26 mm | 1/2 | 75 | blue | | | | DS13PA | | | |
| 26 mm | 1/2 | 90 | blue | | | DS11PA | | | | |
| 30 mm | 1/2 | 75 | blue | | | | | | DS21PA | |
| | | 90 | blue | | | | DS05PA | | DS14PA | |
| 40 mm | 1/2 | 75 | blue | | | | | | DS07PA | DS22PA |
| | | 100 | blue | | | | | DS06PA | | |
| 45 mm | 1/2 | 100 | blue | | | | | DS08PA | | |
| | | 200 | blue | | | | | | DS04PA | |
| 45 mm | 1/2 | 200 | blue | | | | | | DS20PA | |
| 45 mm heavy | 1/2 | 200 | blue | | | | | | | DS23PA |

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting P needle blade

| length | needles curvature | | thread length (cm) color | | | | | | | | | |
|--------|----------------------|--|--------------------------------|-----|-----|--------|--------|--------|---|---|--|--|
| | | | | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | | |
| 13 mm | 3/8 | | 75 blue | | | DM12PA | | | | | | |
| 16 mm | 3/8 | | 45 blue | | | DM10PA | DM01PA | | | | | |
| | | | 45 black | | | DM09PA | | | | | | |
| | | | 75 blue | | | DM11PA | DM02PA | DM03PA | | | | |
| 19 mm | 3/8 | | 45 blue | | | | DM04PA | | | | | |
| | | | 75 blue | | | | DM05PA | DM06PA | | | | |
| 24 mm | 3/8 | | 45 blue | | | | | DM08PA | | | | |
| | | | 75 blue | | | | | DM07PA | | | | |

notes

NYLON (polyamide)

surgical sutures, non-absorbable, monofilament

LIGATURES

| threads length | color thread | 7/0 | 6/0 | 5/0 | 4/0 | 3/0 | 2/0 | 0 | 1 | 2 |
|----------------|--------------|-----|-----|-----|-----|--------|--------|---------|---------|---------|
| 2 x 75 cm | blue | | | | | | P05PA | | P06PA | |
| 2 x 100 cm | blue | | | | | P011PA | P038PA | P012PA | P013PA | P014PA |
| 3 x 45 cm | blue | | | | | P59PA | P49PA | P60PA | | |
| 6 x 45 cm | blue | | | | | P031PA | P032PA | P25PA | P028PA | P026PA |
| 75 cm | blue | | | | | | | | P03PA | P04PA |
| 100 cm | blue | | | | | | P07PA | P08PA | P09PA | P010PA |
| 150 cm | blue | | | | | P016PA | P020PA | P017PA | P018PA | P019PA |
| 150 cm in reel | blue | | | | | | | Pr017PA | Pr018PA | Pr019PA |
| 180 cm | blue | | | | | | P46PA | P45PA | P43PA | P44PA |

notes